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- 1) 山崎 哲, 鈴木典子, 後藤宏美, 他. APTT の現状と標準化に向けた課題. 生物試料分析. 2009 ; **32** (5) : 365–370.
- 2) 向出佳恵. マルチウェーブ検出方式を用いた全自動血液凝固測定装置 CS-2000*i*/2100*i* の特徴. 生物試料分析. 2009 ; **32** (5) : 393–400.
- 3) Toh CH, Giles AR. Waveform analysis of clotting test optical profiles in the diagnosis and management of disseminated intravascular coagulation (DIC). Clin and Lab Haematol. 2002 ; **24** (6) : 321–327.
- 4) Toh CH, Ticknor LO, Downey C, et al. Early identification of sepsis and mortality risks through simple, rapid clot-waveform analysis. Implications of lipoprotein-complexed C reactive protein formation. Intensive Care Med. 2003 ; **29** (1) : 55–61.
- 5) Smith EY, Charles LA, Van Cott EM. Biphasic activated partial thromboplastin time waveform and adverse events in non-intensive care unit patients. Am J Clin Pathol. 2004 ; **121** (1) : 138–141.
- 6) Toh CH, Samis J, Downey C, et al. Biphasic transmittance waveform in the APTT coagulation assay is due to the formation of a Ca⁺⁺-dependent complex of C-reactive protein with very-low-density lipoprotein and is a novel marker of impending disseminated intravascular coagulation. Blood. 2002 ; **100** (7) : 2522–2529.
- 7) Chopin N, Floccard B, Sobas F, et al. Activated partial thromboplastin time waveform analysis : a new tool to detect infection? Crit Care Med. 2006 ; **34** (6) : 1654–1660.
- 8) Ahmed NZ, Simona C, Nuffelen V, et al. Combination of biphasic transmittance waveform with blood procalcitonin levels for diagnosis of sepsis in acutely ill patients. Crit Care Med. 2008 ; **36** (5) : 1507–1512.
- 9) Dempfle CE, Borggreffe M. The hidden sepsis marker : aPTT waveform analysis. Thromb Haemost. 2008 ; **100** (1) : 9–10.
- 10) 前田奈弥, 山本ゆか子, 鈴木敦夫, 他. 全自動血液凝固測定装置 CN-6000 における初期反応異常低減機能の有効性に関する検証. 日本検査血液学会雑誌 検査と血液. 2019 ; 第 20 巻 学術集会号 : S90.
- 11) Brandt JT, Triplett DA, Musgrave K, et al. The sensitivity of different coagulation reagents to the presence of lupus anticoagulants. Arch Pathol Lab Med. 1987 ; **111** (2) : 120–124.
- 12) Denis-Magdelaine A, Flahault A, Verdy E. Sensitivity of sixteen APTT reagents for the presence of lupus anticoagulants. Haemostasis.

- 1995 ; **25** (3) : 98–105.
- 13) Kitchen S, Carwright I, Woods TA, et al. Lipid composition of seven APTT reagents in relation to heparin sensitivity. *Br J Haematol.* 1999 ; **106** (3) : 801–808.
- 14) Lawrie AS, Kitchen S, Purdy G, et al. Assessment of Actin FS and Actin FSL sensitivity to specific clotting factor deficiencies. *Clin Lab Haematol.* 1998 ; **20** (3) : 179–186.

Evaluation of Reduced Error Rate by Improved Early Reaction Error Function for the Automated Blood Coagulation Analyzer CN-6000

Tomoe YAMASHITA^{*1}, Yoko OKAZAKI^{*1}, Takamitsu IMANISHI^{*1}, Yuka TABUCHI^{*2}, Hiroshi KURONO^{*2}, Ryosuke NARISADA^{*3}, Hironori KATSUMI^{*2} and Jun SAEGUSA^{*1}

^{*1} Department of clinical laboratory, Kobe University Hospital, 7-5-2 Kusunoki-cho, Chuo-ku, Kobe, 650-0017, Japan

^{*2} System Engineering, Sysmex Corporation

^{*3} Medical Affairs, Sysmex Corporation

In a measurement of activated partial thromboplastin time (APTT), a sample, which has a biphasic wave showing clear clotting reaction after a primary change in optical turbidity, can seldom be observed. The biphasic wave is expected to be a biomarker. On the Automated Blood Coagulation Analyzers CS-5100, CS-2400, CS-2500, CS-2000i, CS-2100i and CS-1600 (hereinafter referred to as CS series: Sysmex Corporation), the Early Reaction Error (ERE) function was implemented to prevent erroneous reporting of clotting time. As a result, a clotting time was not output when a result was flagged with ERE. Therefore, a reduction in the ERE rate was required.

On the Automated Blood Coagulation Analyzers CN-3000, CN-6000, CN-3500 and CN-6500 (hereinafter referred to as CN series: Sysmex), in addition to a conventional ERE function, a new analysis feature which detects biphasic waves was implemented. These two functions were evaluated in this study for the reduction in the ERE rate.

Analyzer used was CN-6000. Reagents used were Thrombocheck APTT-SLA, RevohemTM APTT SLA and Actin FSL (all from Sysmex). Clinical samples were collected at Kobe University Hospital and those with an ERE flag received during routine APTT measurements were used. Also, plasma samples spiked with CRP and VLDL to initiate biphasic waves were used. Correlation was evaluated with KC4 Delta - semi automated coagulation analyser (DS-MEDICAL inc.) which utilizes a physical detection method.

The new analysis feature showed improved ERE rate to all reagents used. The evaluation results had excellent correlation with KC4 Delta and showed the validity of a clotting time calculated using the new analysis feature. This new feature installed on the CN series is expected to report the clotting time for biphasic waves which was not reportable on the CS series.

Key Words

CN-6000, Optical Detection Method, APTT, Clot Waveform, Biphasic Waveform