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- 1) Bruce L Evatt, et al. : Fundamental Diagnostic Hematology, The Bleeding and Clotting Disorders (second edition). World Health Organization, 1992.
- 2) NCCLS. Determination of factor coagulant activities; approved guideline H48-A. Wayne, PA: 1997.
- 3) NCCLS: Procedure for the determination of fibrinogen in plasma; approved guideline H30-A2 second edition. Wayne, PA: 2001.
- 4) Donna M Corriveau, et al. : Hemostasis and Thrombosis in the Clinical Laboratory, J B Lippincott Company, USA, 1988.
- 5) K N Williams, et al. : A Computer Program for the Analysis of Parallel-Line Bioassay of Clotting Factors. British Journal of Haematology, 31: 13, 1975.
- 6) T B L Kirkwood, et al. : Biometric principles in clotting and clot lysis assays. Clin lab Haemat 2: 155-167, 1980.

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Severe = less than 1% activity
Moderately severe = 1% - 5% activity
Mild severe = 5% - 25% activity

Evaluation of Factor VIII Assay using Multi Dilution Assay (MDA) on CA Analyzers

Nobuo ARAI, Yuko OGURA, Toshikazu KASAI, and Tsuyoshi ISHII

Scientific Division, Sysmex Corporation
4-4-4 Takatsukadai, Nishi-ku, Kobe 651-2271.

SUMMARY

The Multi Dilution Assay (MDA) is a function to analysis the Factor VIII Activity, and to check the presence of inhibitor or severe factor deficiency. MDA Line is the straight line which connects the plotted measured values by auto doubling dilution of a patient sample. Comparing the shape, linearity, slope and parallel with the calibration curve, the inhibitor or factor deficiency can be suspected.

In current system software, the Slope Ratio (SR) has been added as a new function. This index is to show the degree of parallelism of the MDA Line with the calibration curve. We investigated the usefulness of SR using Factor VIII Inhibitor Plasma and Factor VIII Deficient Plasma and reported on the results as follows.

Key Words Multi Dilution Assay (MDA), Factor VIII Assay, APTT, Slope Ratio (SR)