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- 1) Freidwald WT, Levy RI, Fredrickson DS : Estimation of the Concentration of Low-Density Lipoprotein Cholesterol in Plasma, Without Use of the Preparative Ultracentrifuge. Clin Chem, 18 : 499 ~ 502, 1979.
- 2) Havel RJ, Eder HA, Brangdon JH : The distribution and chemical composition of ultracentrifugally separated lipoproteins in human serum. J Clin Invest, 34 : 1345 ~ 1353, 1955.
- 3) Demacker PN, et al.: Five methods for determining low-density lipoprotein cholesterol compared. Clin Chem, 30 : 1797 ~ 1800, 1984.
- 4) McNamara JR, et al.: Calculated values for low-density lipoprotein cholesterol in the Assessment of lipid abnormalities and coronary disease risk. Clin Chem, 36 : 36 ~ 42, 1990.
- 5) Warnick GR, et al.: Estimating low-density lipoprotein cholesterol by the Friedwald equation is adequate for classifying patients on the basis of nationally recommended cutpoints. Clin Chem, 36 : 15 ~ 19, 1990.
- 6) 岸浩司, 他 : LDL-C試薬・KL「コクサイ」の基礎的検討. Sysmex J, 25 : 14 ~ 21, 2003.
- 7) Becher JD, et al.: Measurement of Low-density lipoprotein cholesterol Concentration, In : Rifai N and Warnick GR, ed. Laboratory Measurement of Lipids, Lipoproteins and Apolipoproteins. Washington, DC : American Association of Clinical Chemistry, 107 ~ 123, 1994.
- 8) 日本動脈硬化学会. 動脈硬化性疾患診療ガイドライン2002年版.
- 9) 厚生労働省 : 標準的な検診・保健指導プログラム(暫定版)

Basic Evaluation of LDL-C Reagent • KL

Ikumi MIYAMOTO, Hirokazu KAKUDA, Masako KADOTA, Nobuo ARAI,
Hirofumi HAGINO, Kunihiro FUNAKOSHI and Keiji FUJIMOTO

Scientific Affairs, Sysmex Corporation,
1-3-2, Murotani, Nishi-ku, Kobe 651-2241

SUMMARY

Sysmex Corporation developed LDL-C Reagent•KL which was a direct method assay kit for LDL cholesterol measurement in 2002. This time, we evaluated the new kit of which shelf life had improved.

In the evaluation of within-run reproducibility, we used the control sera of 3 concentration. CV_s were shown in 0.63~1.42%. The linearity was observed up to 400mg/dL. Hemoglobin (500mg/dL), bilirubin (30mg/dL) and ascorbic acid (100mg/dL) did not interfere in the assay. The correlation coefficient using fresh sera (n=53) between the present kit and the improved kit was 0.998 and the regression formula was $y=0.997x+1.33$. On the other hand, the correlation coefficient between the improved kit and -quantification method (BQ method) was 0.988 and the regression formula was $y=1.005x+1.07$. It was confirmed that the shelf life of kit was extended to 12 months from 8 months.

Key Words LDL cholesterol, BQ method