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- 1) James GP, Harrison R. Creatine kinase isozyme of mitochondrial origin in human serum. *Clin Chem.* 1979; 25(6):943-947.
- 2) Payne RM, Strauss AW. Expression of the mitochondrial creatine kinase genes. *Mol Cell Biochem.* 1994; 133-134: 235-243.
- 3) Payne RM, Strauss AW. Developmental expression of sarcomeric and ubiquitous. *Mol Cell Biochem.* 1994; 1219(1):33-38.
- 4) Anfous K, et al. Mitochondrial creatine kinase isoform expression does not correlate with its mode of action. *Biochem J.* 1997; 322:73-78.
- 5) Haas RC, Strauss AW. Separate Nuclear Genes Encode Sarcomere-specific and Ubiquitous Human Mitochondrial Creatine Kinase Isoenzymes. *J Biol Chem.* 1990; 265: 6921-6927.
- 6) 星野 忠. ミトコンドリアCKの測定法. *検査と技術.* 2000; 28(13):1499-1504.
- 7) Kanemitsu F, et al. Characterization of two types of mitochondrial creatine kinase isolated from normal human

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Basic Performance of "L-System CK", Reagents for Creatine Kinase Measurement and "L-System CK-MB", Reagents for Creatine Kinase Isozyme Measurement

Ikumi Miyamoto, Hirokazu Kakuda, Nobuo Arai, Hirofumi Hagino and Kunihiro Funakoshi

Scientific Affairs, Sysmex Corporation, 1-3-2, Murotani, Nishi-ku, Kobe 651-2241

SUMMARY

Creatine kinase (CK) activity is widely used as a diagnostic parameter of muscle disorder and cardiac disease. CK has two major subunits of M and B, and there are three major isoenzymes: CK-BB (as found in the brain), CK-MM (as found in the skeletal muscle), and CK-MB (in the cardiac muscle). CK-MB is found mainly in myocardial tissue, so this isoenzyme is useful in a diagnosis of myocardial infarction. Recently, we developed new reagents for CK and CK-MB activity measurement avoiding the interference of hemoglobin and mitochondrial CK; L-System CK and L-System CK-MB.

In this paper, we report the basic performance of both new reagents.

L-System CK and L-System CK-MB showed the excellent results in within-run reproducibility and linearity. Moreover, there is no influence of hemoglobin up to 500mg/dL on the measurement results of both reagents. Anti-mitochondria CK antibody which was contained in L-System CK-MB was able to inhibit 98.3% of sarcomeric mitochondrial CK and 96.3% of ubiquitous mitochondrial CK. As for method comparison, the correlation coefficient and regression equation between L-System CK-MB and the conventional reagent, CK-MB reagent · L "kokusai" were 0.900 and $y = 0.799x - 8.08$ respectively and the discrepancy between measurement result of the two reagents was found in some samples. However, on the correlation between L-System CK-MB and electrophoresis assay, a good result was shown ($r = 0.971$, $y = 0.909x + 5.1$). ROC analysis in cardiac disease showed a larger AUC (Area Under the Curve) of L-System CK-MB (0.906) than the conventional reagent (0.802).

From these observations, it was shown that L-System CK and L-System CK-MB have better performance, and L-System CK-MB is able to measure CK-MB activity more specifically than the conventional reagent.

Key Words CK, CK-MB, Mitochondrial CK (MtCK), Anti-human MtCK antibody, Immunological inhibition method