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
※1 HGMD® is a registered trademark of Cardiff University

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1. 1,000ゲノム
2. gnomAD
3. dbSNP
4. GEM-J WGA

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31 Jul 2023 v1
posi_con

PrismGuide IRD Panel Summary Report (Prioritized Variants)

- ① **QC Status: Pass**
- ② Order ID: posi_con
- ③ **Sample Information**

Clinical Sample ID	posi_con
--------------------	----------
- ④ **SNV/Indel Information**

#	Gene Name	Inheritance Pattern	Mutant Allele Frequency		Mutation Content		Consequence	
			Frequency (Read)	Zygosity	HGVSc	HGVSp	ACMG Classification	ACMG labels
1	ABCA4	Recessive	53.29% (721/1353)	Heterozygous	c.5318C>T	p.Ala1773Val	Pathogenic	PS1,PS3,PM2,PP2,PP3
2	ABCA4	Recessive	53.70% (610/1136)	Heterozygous	c.1760+2T>G	-	Pathogenic	PVS1,PS3,PM2,PP3,PP5
3	RPE65	Recessive, Dominant	53.73% (663/1234)	Heterozygous	c.1154C>T	p.Thr385Met	Pathogenic	PS1,PS3,PM2,PP2,PP3,BP4
4	RPGR	X-linked	98.68% (149/151)	Hemizygous	c.2838_2839del	p.Glu947GlyfsTer131	Pathogenic	PVS1,PS3,PM2
5	CEP290	Recessive	24.35% (158/649)	Heterozygous	c.2991+1655A>G	-	Likely pathogenic	PS3,PM2
- ⑤ **SNV/Indel Report:**
 1. ABCA4のAla1773Valはヘテロであり、潜在(劣性)の病的変異である。
 2. ABCA4のc.1760+2T>Gはヘテロであり、潜在(劣性)の病的変異である。
 3. RPE65のThr385Metはヘテロであり、潜在(劣性)または顕性(優性)の病的変異である。
 4. RPGRのGlu947GlyfsTer131はホモであり、X-linkedの病的変異である。
 5. CEP290のc.2991+1655A>Gはヘテロであり、潜在(劣性)の病的変異の可能性がある。
- ⑥ SNV/Indel Notes:
- ⑦ Sample Notes:
- ⑧ Report Date: Confirmation Sign:

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31 Jul 2023 v1	posi_con
PrismGuide IRD Panel Sequencing Report	
①	Order ID: posi_con
②	Sample Information
	Clinical Sample ID: posi_con
	Customer-provided Sample Sex: Male
③	Sequencing Analysis Information
	Panel: PrismGuide IRD Panel
	Reagent: PrismGuide IRD Panel reagent kit
	Analysis Execution Date: 31 Jul 2023 14:25:25
	On Target Reads (%): 73.86
	On And Near Target Reads (%): 82.55
	Off Target Reads (%): 17.45
	Mapped Reads (%): 99.37
	Duplicated Reads (%): 6.33
	Reads Mapping Quality 0 (%): 3.53
	Average Quality: 36.5
	Average Insert Size: 195.1
	Insert size std: 55.3
	Mean Target Coverage: 588.33
	Targets Not Covered: 4
	Evenness: 83.83
	Uniformity: 1.573
	Derivative Log Ratio Score for PrismGuide IRD Panel System: 0.133
	Pipeline-calculated Sample Sex: Male
	Low Coverage (>=100x): 98.89
	Medium Coverage: 65.27
	High Coverage: 2.14
④	Data Analysis Information
	Analysis Protocol: IRD Protocol v1
	Pipeline Version: IRD Pipeline v2.4.3
	Enable CNV Calling: Yes
	Import SNV Results: Yes
	CNV Calling Mode: Yes
	Flanking Region (bp): 10
	Allele Balance Priors Off: Yes
	Minimum Alt Count: 5
	Minimum Base Quality: 20
	Minimum Total Read Count: 20
	Minimum Mapping Quality: 30
	Max Read Mismatch Fraction: 0.04
	Threshold for Targets Not Covered: 1
	Germline Minimum Alt Fraction: 0.2
	Germline Segmentation p-value Threshold: 0.01
	Genetic aberration selection criteria (SNV, InDel): Default SNV Filter
	Genetic aberration selection criteria (CNV): Default CNV Filter
	Analysis Program Version: 0.3.20

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6 CNV Information (as Reference)

1	Variant (ISCN)	seq[GRCh38] 1q31.3(195134497_197328497)x1
	Variant Type	Deletion
	Location	1:195134497-197328497
	Copy Number	1
	Log Ratio	-0.617
	Genes	CRB1 (exons 1-2/12)
	User Classification	Unclassified
2	Variant (ISCN)	seq[GRCh38] 2q13q13(111996876_111997026)x1
	Variant Type	Deletion
	Location	2:111996876-111997026
	Copy Number	1
	Log Ratio	-1.055
	Genes	MERTK (intronic)
	User Classification	Unclassified
3	Variant (ISCN)	seq[GRCh38] 4p15.32p15.32(15979817_15979967)x1
	Variant Type	Deletion
	Location	4:15979817-15979967
	Copy Number	1
	Log Ratio	-1.425
	Genes	PROM1 (exon 25/28)
	User Classification	Unclassified
4	Variant (ISCN)	seq[GRCh38] 5q14.3q14.3(90724940_90725233)x1
	Variant Type	Deletion
	Location	5:90724940-90725233

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7 Appendix 1: Disease Category Referred From RetNet

Gene	Associated Diseases
ABCA4	Cone or cone-rod dystrophy, autosomal recessive; Macular degeneration, autosomal recessive; Retinitis pigmentosa, autosomal recessive
ADGRV1	Usher syndrome, autosomal recessive
AIP1	Cone or cone-rod dystrophy, autosomal dominant; Leber congenital amaurosis, autosomal recessive
BEST1	Macular degeneration, autosomal dominant; Retinitis pigmentosa, autosomal dominant; Retinitis pigmentosa, autosomal recessive; Other retinopathy, autosomal dominant; Other retinopathy, autosomal recessive
C8orf37	Bardet-Biedl syndrome, autosomal recessive; Cone or cone-rod dystrophy, autosomal recessive; Retinitis pigmentosa, autosomal recessive
CA4	Retinitis pigmentosa, autosomal dominant
CACNA1F	Cone or cone-rod dystrophy, X-linked; Congenital stationary night blindness, X-linked; Other retinopathy, X-linked
CDH23	Deafness alone or syndromic, autosomal recessive; Usher syndrome, autosomal recessive
CDHR1	Cone or cone-rod dystrophy, autosomal recessive
CEP290	Bardet-Biedl syndrome, autosomal recessive; Leber congenital amaurosis, autosomal recessive; Syndromic/systemic diseases with retinopathy, autosomal recessive
CERKL	Cone or cone-rod dystrophy, autosomal recessive; Retinitis pigmentosa, autosomal recessive
CFAP410	Cone or cone-rod dystrophy, autosomal recessive
CHM	Other retinopathy, X-linked
CLRN1	Retinitis pigmentosa, autosomal recessive; Usher syndrome, autosomal recessive
CNGA1	Retinitis pigmentosa, autosomal recessive
CNGA3	Cone or cone-rod dystrophy, autosomal recessive; Other retinopathy, autosomal recessive
CNGB1	Retinitis pigmentosa, autosomal recessive
CNGB3	Cone or cone-rod dystrophy, autosomal recessive; Other retinopathy, autosomal recessive
CRB1	Leber congenital amaurosis, autosomal recessive; Retinitis pigmentosa, autosomal recessive; Other retinopathy, autosomal dominant
CRX	Cone or cone-rod dystrophy, autosomal dominant; Leber congenital amaurosis, autosomal dominant; Leber congenital amaurosis, autosomal recessive; Retinitis pigmentosa, autosomal dominant
CYP4V2	Retinitis pigmentosa, autosomal recessive; Other retinopathy, autosomal recessive
DHDDS	Retinitis pigmentosa, autosomal recessive
DRAM2	Macular degeneration, autosomal recessive

8 Appendix 2: Criteria for classifying pathogenic and benign variants based on ACMG guidelines

Very strong evidence of pathogenicity	
PVS1	Null variant (nonsense, frameshift, canonical ± 1 or 2 splice sites, initiation codon) in a gene where loss of function (LOF) is a known mechanism of disease. SOURCE=VEP Impact
Strong evidence of pathogenicity	
PS1	Same amino acid change as a previously established pathogenic variant regardless of nucleotide change. SOURCE=ClinVar, HGMD, Custom Variant List
PS3	Well-established in vitro or in vivo functional studies supportive of a damaging effect on the gene or gene product. SOURCE=ClinVar, HGMD
Moderate evidence of pathogenicity	
PM1	Located in a mutational hot spot and/or critical and well-established functional domain without benign variation. SOURCE=ClinVar, PFAM
PM2	Absent from controls if dominant or with a frequency $\leq 0.5\%$ if recessive in gnomAD database. SOURCE=gnomAD Allele Frequency
PM4	Protein length changes due to in-frame deletions/insertions in a non-repeat region or stop-loss variants. SOURCE=VEP consequence, Repeat Masker
PM5	Novel missense change at an amino acid residue where a different missense change determined to be pathogenic and has been seen before. SOURCE=ClinVar, HGMD, Custom Variant List
Supporting evidence of pathogenicity	
PP2	Missense variant in a gene that has a low rate of benign missense variation and where missense variants are a common mechanism of disease. SOURCE=ClinVar
PP3	Multiple lines of computational evidence support a deleterious effect on the gene or gene product (conservation, evolutionary, splicing impact, etc.). SOURCE=SIFT, BLOSUM62, dbcsSNV, GERP
PP5	Reputable source recently reports variant as pathogenic but the evidence is not available to the laboratory to perform an independent evaluation. SOURCE=ClinVar, HGMD, Pathogenic Variants from PMID:31213501
Stand-Alone evidence of benign impact	
BA1	Allele frequency is above 5% in gnomAD database. SOURCE=gnomAD Allele Frequency
Strong evidence of benign impact	
BS1	Allele frequency is greater than 0.5% for disorder. SOURCE=gnomAD Allele Frequency
BS3	Well-established in vitro or in vivo functional studies show no damaging effect on protein function or splicing. SOURCE=ClinVar
Supporting evidence of benign impact	
BP1	Missense variant in a gene for which primarily truncating variants are known to cause disease. SOURCE=ClinVar
BP3	In-frame deletions/insertions in a repetitive region without a known function. SOURCE=VEP consequence, Repeat Masker, PFAM
BP4	Multiple lines of computational evidence suggest no impact on gene or gene product (conservation, evolutionary, splicing impact, etc.). SOURCE=SIFT, BLOSUM62, dbcsSNV, GERP
BP6	Reputable source recently reports variant as benign but the evidence is not available to the laboratory to perform an independent evaluation. SOURCE=ClinVar, Custom Variant List, Benign Variants from PMID:31213501

9 Appendix 3: Rules for combining criteria to classify variants based on ACMG guidelines

- Pathogenic
 1. 1 Very Strong (PVS1) AND
 1. ≥ 1 Strong (PS1,PS3) OR
 2. ≥ 2 Moderate (PM1,PM2,PM4,PM5) OR
 3. 1 Moderate (PM1,PM2,PM4,PM5) and 1 Supporting (PP2,PP3,PP5) OR
 4. ≥ 2 Supporting (PP2,PP3,PP5)
 2. ≥ 2 Strong (PS1,PS3) OR
 3. 1 Strong (PS1,PS3) AND
 1. ≥ 3 Moderate (PM1,PM2,PM4,PM5) OR
 2. 2 Moderate (PM1,PM2,PM4,PM5) AND ≥ 2 Supporting (PP2,PP3,PP5) OR
 3. 1 Moderate (PM1,PM2,PM4,PM5) AND ≥ 3 Supporting (PP2,PP3,PP5)
- Likely Pathogenic
 1. 1 Very Strong (PVS1) AND 1 Moderate (PM1,PM2,PM4,PM5) OR
 2. 1 Strong (PS1,PS3) AND 1-2 Moderate (PM1,PM2,PM4,PM5) OR
 3. 1 Strong (PS1,PS3) AND ≥ 2 (PP2,PP3,PP5) OR
 4. ≥ 3 Moderate (PM1,PM2,PM4,PM5) OR
 5. 2 Moderate (PM1,PM2,PM4,PM5) AND ≥ 2 Supporting (PP2,PP3,PP5) OR
 6. 1 Moderate (PM1,PM2,PM4,PM5) and ≥ 3 Supporting (PP2,PP3,PP5)
- Benign
 1. 1 Stand-Alone (BA1) OR
 2. ≥ 2 Strong (BS2,BS3)
- Likely Benign
 1. 1 Strong (BS1, BS3) AND 1 Supporting (BP1,BP3,BP4,BP6,BP7) OR
 2. ≥ 2 Supporting (BP1,BP3,BP4,BP6,BP7)
- Uncertain Significance
 1. Variants that do not meet the above criteria
 2. Conflicting classification labels (eg Likely Pathogenic, Likely Benign) which are listed

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12 Appendix 6: QC Metrics

Name	Definition
Mapped Reads (%)	The percentage of total reads that are aligned to the reference genome
Duplicated Reads (%)	The percentage of total reads that are marked duplicated
Mean Target Coverage	The mean target coverage is defined as the sum of base counts (coverage) within the targeted regions of the panel divided by the total length of these regions. All base counts are calculated from filtered, mapped and de-duplicated reads
Low Coverage (>=100x)	The percentage of bases that have a coverage above or equal to 100
Medium Coverage	The percentage of bases that have a coverage above or equal to 500
High Coverage	The percentage of bases that have a coverage above or equal to 1000
Evenness	The evenness score represents the fraction of whole-sequencing throughput that is correctly distributed over the targeted regions of the panel
Average Quality	The average base Phred quality score of all reads
Average Insert Size	The average absolute template length for paired and mapped reads

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Name	Pass
Mapped Reads (%)*	>=95.0
Mean Target Coverage	>=240.0
Low Coverage (>=100x)	>=95.0
Evenness	>=75.0
Targets Not Covered*	<=30.0
On Target Reads (%)	>=60.0

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